

DIVATOR MKII SET

Certificate no: DSC-2026-000608

Service date: 06.04.2026

Valid until: 06.04.2027

Result: **Approved**

Equipment details

Customer: Sjøentreprenøren AS
Producer: Interspiro
Serial no: 10380691
Service interval: 12 months

Location: Alu kontainer
Model: Divator MKII
Build year: -

NOTE

Tested with Labtec PRO

This is a certificate confirmation, the documentation has a separate file.

Technician

Serviced by: Daniel Kristensen

Signature



test report

test report: 455
 device no.: AL001
 type of Yearly test
 test system: PRO 3102.24.325

workshop: DS BERGEN AS		person:	Daniel Kristensen
street:	Storebotn 55G	phone:	91649162
place:	5309 Kleppestø	fax:	
		email:	Daniel@dsbergen.no
customer: Sjøentreprenøren AS		person:	Ingmar Aune
street:		phone:	
place:	Bergen	fax:	
		email:	ingmar.aune@ncc.no
device data:	deadlines:	device-info:	service:
device no: AL001	1 y test 06.04.2027 1 y srv 06.04.2027 6 y srv 14.05.2031		- 1 year service - cleaning and disinfection
type: Interspiro Divator NP			
class: mask NP			
location: Alu kontainer			

Divator mask, m/ Kirby blokk.
 Uoriginal BO.

tests/measurements	unit	min	max	actual	result
visual test (0 = ok) input mask: (0 = classic, > 0 = extended/list...)			0	0	OK
mask fitted approx. 6 sec head expansion Mask volume 0,2 l				0	OK
leakage test w/ open exhalation valve 0,25 min at p=-10 mbar	mbar		7,5	2,5	OK

test executed: 06.04.2026 13:30
 by: Daniel Kristensen
 test signed: 06.04.2026 13:36
 by: Daniel Kristensen

result:
device ready for use

LabNet V: 2025.12.2.0

next test on: 06.04.2027

workshop: DS BERGEN AS		person:	Daniel Kristensen
street: Storebotn 55G		phone:	91649162
place: 5309 Kleppestø		fax:	
		email:	Daniel@dsbergen.no
customer: Sjøentreprenøren AS		person:	Ingmar Aune
street:		phone:	
place: Bergen		fax:	
		email:	ingmar.aune@ncc.no
device data:	deadlines:	device-info:	service:
device no: 1038069L	1 y test 06.04.2027 1 y srv 06.04.2027 6 y srv 14.05.2031		1 year service cleaning and disinfection
type: Interspiro BV Divator NP			
class: demand valve NP			
location: Alu kontainer			

Divator mask, m/ Kirby blokk.
Uoriginal BO.

tests/measurements	unit	min	max	actual	result
visual test (0 = ok) input mask: (0 = classic, > 0 = extended/list...)			0	0	OK
demand valve opening pressure at 10 l/min suction	mbar	-4	-0,5	-2,4	OK
inhalation resistance at 40 x 2,5 l/min, 5 hubs HP .>. bar	mbar	-8	-0,1	-3,9	OK
exhalation valve opening pressure at 10 l/min suction	mbar	2	5	2,4	OK
safety valve opening pressure	bar	12	20	15,5	OK
pressure release of test equipment (HP, MP and LP)				0	OK

test executed: 06.04.2026 13:30
 by: Daniel Kristensen
 test signed: 06.04.2026 13:36
 by: Daniel Kristensen

result:
device ready for use

LabNet V: 2025.12.2.0

next test on: 06.04.2027